

50 pcs

non-aseptic non-reusability

Epidemic emergency product, It is only used for general protection of non-medical personnel during the epidemic prevention and control in this region

non-woven fabrics,ear loop 17.5cm×9.5cm 3layers

# DISPOSABLE FACE MASK

50 pcs

# DISPOSABLE FACE MASK

non-woven fabrics,ear loop 17.5cm×9.5cm 3layers

non-reusability



**[Product name]** Disposable Face Mask **50PCS/BOX**  
**[Main material]** Non-woven and melt blown  
**[Specification type]** 17.5\*9.5cm,3 layers mask,elastic rope for ear; metal strips for nose wire  
**[Product protection level]** II

**[Instruction for using]**  
1.Nose clip up, hold one ear curtain with both hands  
2.Put on the mask, put the chin in the mask, buckle the ear bag behind the ear with both hands, and make the mask stick to the face  
3.Use the forefinger and middle finger of both hands to press and adjust the nasal clamp to stick to the bridge of the nose

**[Be attention]** If dizziness, nausea and dyspnea occur during use, please replace the product. This product is disposable and will be destroyed after use.



# DISPOSABLE FACE MASK

non-woven fabrics,ear loop 17.5cm×9.5cm 3layers

non-reusability



50 pcs

**Lot:** 2020 Mar 1st  
**Expired date:** 2025 Mar 1st  
**Package:** 50pcs  
**Distributor:** OUDA INTERNATIONAL LIMITED  
**Factory:** Ganzhou Zhongke TUOYOU DA Intelligent Equipment Technology Co.,Ltd  
**Address:** Room A01,10th Floor,JinZhou international business center,No.899,Jiefang North Road,Yuexiu District,Guangzhou  
City

MADE IN P.R.C